

Silicone Drains and Wound Drains



Silicone Drains

Bellovac®
Exudrain®

The perfect fit
in every way



Bellovac®



Large Volume Drainage

Average initial negative pressure when activated: 115 mmHg/15 kPa
Volume bellows: 220 ml
Volume bag: 700 ml

Surgical Procedures:

Orthopaedics
Gynaecology
General Surgery

Ref No Bx/Pk Description


68390	x15	Bellovac® Set FG 10
68391	x15	Bellovac® Set FG 14
68392	x15	Bellovac® Set FG 18
68393	x15	Bellovac® Set FG 10 with trocar
68394	x15	Bellovac® Set FG 14 with trocar
68395	x15	Bellovac® Set FG 18 with trocar

Y-Connectors

Ref No Bx/Pk Description

68376	x25	Y-Connectors w/ FG10 catheter
68377	x25	Y-Connectors w/ FG14 catheter
68378	x25	Y-Connectors w/ FG18 catheter
68385	x25	Y-Connectors w/ FG10 catheter and trocar
68386	x25	Y-Connectors w/ FG14 catheter and trocar
68387	x25	Y-Connectors w/ FG18 catheter and trocar

Exudrain®



Medium Volume Drainage

Average initial negative pressure when activated: 75 mmHg/10 kPa
Volume bulbs: 110 ml
Volume bag: 700 ml

Surgical Procedures:

General Surgery
Breast Surgery, e.g. mastectomies
Wound drain in situ

Ref No Bx/Pk Description

68406	x10	Exudrain® Set FG 10
68407	x10	Exudrain® Set FG 14
68408	x10	Exudrain® Set FG 18
68409	x10	Exudrain® Set FG 10 with trocar
68410	x10	Exudrain® Set FG 14 with trocar
68411	x10	Exudrain® Set FG 18 with trocar
68419	x10	Exudrain® Mini FG 8 with trocar
68412	x10	Exudrain® Mini FG 8

Exchange Bags

Ref No Bx/Pk Description

68337	x50	Exchange Bag
68480	x25	Sterile Exchange Bag

Silicone Range



Silicone Drainage

Flat Fluted Silicone, Round Fluted Silicone, Flat Silicone and Round Silicone Catheter range.
For safe easy drainage use with Bellovac® and Exudrain® drainage systems.

Surgical Procedures:

Plastic Surgery, ENT Surgery
General Surgery, Neuro Surgery

Ref No Bx/Pk Description

68463	x15	Fluted Flat Silicone Drain 7mm/FG14 with Trocar
68464	x15	Fluted Flat Silicone Drain 10mm/FG14 with Trocar
68460	x15	Fluted Silicone Drain FG10 with Trocar
68461	x15	Fluted Silicone Drain FG14 with Trocar
68462	x15	Fluted Silicone Drain FG18 with Trocar
68361	x15	Silicone Catheter (Flat) with Trocar 14FG/7mm
68362	x15	Silicone Catheter (Flat) with Trocar 14FG/10mm
68367	x15	Silicone Catheter (Round) with Trocar 10FG
68368	x15	Silicone Catheter (Round) with Trocar 14FG
68369	x15	Silicone Catheter (Round) with Trocar 18FG