

# Bamboo

MEDIPLAST

Stretchable tubular  
bandage, made of bamboo.



LONG-TERM FRESHNESS



WONDERFULLY SOFT




EXCELLENT HEAT REGULATION



ECO-FRIENDLY





## Danafast DS

MEDIPLAST offers DANAFAST Double Stretch, a tubular elastic bandage that stretches in length and width and is thus easy to apply. Can be cut in desired length and can be applied in single or double layer, no additional fixation is required. Since the product is knitted it provides a good air-permeability. Colour coded with a woven line which makes choice of size easy.



## Danagrip

MEDIPLAST offers DANAGRIP, an elastic support and fixation bandage with a light compression effect. Holes and gashes can be cut out. Can be used in single or double layers depending on desired compression. Since the product is knitted it provides a good air-permeability. Good elasticity. The bandages are available in three different materials, unbleached cotton, bleached cotton, and bamboo.



## Danatube

MEDIPLAST offers DANATUBE, an elastic tubular bandage. Since the product is knitted it provides a good air-permeability. The bandages are available in four different materials, unbleached cotton, bleached cotton, viscose and bamboo.